(EN) (DE) (FR) (ES) (PT) (AR)

ACCU-CHEK® Active

03146324

(EN) **Control**

| IN VITRO DIAGNOSTIC MEDICAL DEVICE – DO NOT INGEST! |
|---|
| Intended Use |
| The control solution is intended for performing control tests on the dedicated test strips and blood glucose meters. It is indicated for self-testing by people with diabetes and for near-patient testing by healthcare professionals. |
| Dedicated test strips and devices are the Accu-Chek Active test strips and the Accu-Chek Active blood glucose meter. |
| Before You Get Started |

Read this package insert, the package insert of the Accu-Chek Active test strips, and the User's Manual of the Accu-Chek Active blood glucose meter before using the control solutions. The User's Manual contains all the information you need to perform a control test and to understand the control results. You can find information on the performance characteristics of the Accu-Chek Active system in the package insert of the Accu-Chek Active test strips.

If you have any questions, contact customer support.

The package insert contains warnings:

A **WARNING** indicates a foreseeable serious hazard.

| ⚠ WARNING |
|--|
| Risk of suffocation |
| This product contains small parts that can be swallowed. Keep the small parts away from small children and people who might swallow small parts. |

Do not ingest or inject the control solution or use the control solution as eye drops.

EUH208 – Contains 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1). May produce an allergic reaction.*

* Regulation (EC) No 1272/2008 (CLP)

The control results must be within the control range shown by the concentration table on the label of the test strip container. The control range specifies the acceptable results for control tests.

If the control result is not within the acceptable control range, it is no longer certain that the meter and test strips are working properly. Do not use the meter and test strips. Contact customer support.

| Contents of the Pack |
|---|
| <ul style="list-style-type: none">1 bottle with control solution with a low glucose concentration (Control 1, blue cap and blue text on the label) 1 bottle with control solution with a high glucose concentration (Control 2, red cap and red text on the label) 1 package insert |

| Additional Materials Required for Performing a Control Test |
|---|
| <ul style="list-style-type: none">Accu-Chek Active blood glucose meter with User's Manual Accu-Chek Active test strips with package insert A clean, dry paper towel |

| Storing and Using Control Solutions Properly |
|---|
| Control solutions, which are not properly stored or used, can deliver results that are outside the specified control range. |
| <ul style="list-style-type: none">The bottles are filled with control solution under low-germ conditions. To avoid contamination, do not touch the tip of the bottle with your fingers or with the test strip. Contaminated control solutions can deliver incorrect control results. Use only control solutions, which are within the use by date. The use by date is printed next to the ☞ symbol on the bottle label and on the cardboard box. Store the control solutions at temperatures between +2 and +30 °C. Control solutions that are stored outside this temperature range may become unusable before they reach their use by date. If you store and handle control solutions as described, the control solutions are valid for 3 months after they were opened for the first time. However, do not use control solutions after the use by date shown next to the ☞ symbol on the bottle. Control solutions must be at ambient temperature when you perform a control test. If you keep the control solutions in a refrigerator, leave them to stand at ambient temperature until they have warmed up before you perform a control test. |

| Ingredients |
|--|
| Composition in percent by weight |
| Control 1: glucose (0.04 %), preservative (0.24 %) in phosphate buffer |
| Control 2: glucose (0.12 %), preservative (0.24 %) in phosphate buffer |

| Discarding the Bottles |
|--|
| All components of the pack can be discarded in domestic waste. |

Discard empty or partly used bottles according to local regulations. For information on how to discard the bottles correctly, contact your local council or authority.

| Reporting of Serious Incidents |
|--|
| For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority. |
| Last Update |
| 2021-05 |

| Customer Support |
|--|
| United Kingdom |
| Roche Diabetes Care Limited |
| Charles Avenue, Burgess Hill |
| West Sussex, RH15 9RY, United Kingdom |
| Accu-Chek Customer Careline ¹⁾ |
| UK Freephone number: 0800 701 000 |
| ROI Freephone number: 1 800 709 600 |
| ¹⁾ calls may be recorded for training purposes |
| Some mobile operators may charge for calls to these numbers. |
| www.accu-chek.co.uk |
| www.accu-chek.ie |

| | |
|------------------|---|
| ⓘ | Consult instructions for use or consult electronic instructions for use |
| ⚠ | Caution, refer to safety-related notes in the package insert accompanying this product. |
| 🌡 | Temperature limit |
| 🕒 | Use by |
| 🕒 | Use-by period after being opened: 3 months |
| 🗑 | All components of the pack can be discarded in domestic waste. Discard empty or partly used bottles according to local regulations. |
| 🏭 | Date of manufacture |
| IVD | In vitro diagnostic medical device |
| 🧑 | Device for self-testing |
| 🧑 | Device for near-patient testing |
| 🏭 | Manufacturer |
| UDI | Unique device identifier |
| REF | Catalogue number |
| LOT | Batch code |
| CE | Complies with the provisions of the applicable EU Legislation |

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