# Implantation set

## Instructions for use

EN DE FR IT NL

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












### (EN) Accu-Chek DiaPort Implantation set

#### Intended Use

The Accu-Chek DiaPort system is intended for continuous intraperitoneal insulin infusion (CIPI). For this purpose, it shall be used together with an external Accu-Chek Spirit Combo insulin pump. The Accu-Chek DiaPort system delivers insulin from the Accu-Chek insulin pump through the infusion set and the catheter into the peritoneal cavity.

**WARNING**  
**Risk of hyperglycaemia**  
The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type 1 or type II diabetes who meet the indication for CIPI therapy. The Accu-Chek DiaPort system is only for use by patients and healthcare professionals specifically trained in the use of the individual components related to their areas of expertise.

**WARNING**  
**Risk of inappropriate therapy**  
The Accu-Chek DiaPort system is designed for intended users only. Do not use the system with patients that do not meet the indications, or where one or more contraindications apply.

The Accu-Chek DiaPort implantation set is a surgical tool kit for the implantation of the Accu-Chek DiaPort system. The product is intended only for the use by surgeons trained in implanting the Accu-Chek DiaPort system, using the tools provided with the set. The implantation must be performed in an operating theater under sterile conditions.

**WARNING**  
**Risk of different types of harm**  
Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explanation of the system.

**WARNING**  
**Risk of inappropriate therapy**  
The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

**WARNING**  
**Risk of infection**  
A damaged packaging may no longer be sterile or may contain a damaged component. Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

**WARNING**  
**Risk of injury**  
If you modify the system or if you use tools that were not provided with the set, you may injure internal organs and tissues.

**WARNING**  
**Risk of hyperglycaemia and infection**  
The use of sterile components after their Use by date can lead to infections

and to abscesses. Do not use components that have passed their Use by date. Before You Get Started Discuss the implantation site with the patient prior to the surgery. Make sure to select a 2nd placement option in case there is a reason as to why the first choice is not possible. Do not use this product for any medication other than U100 short-acting regular human insulin or sterile physiological saline solution. Read the Accu-Chek DiaPort Handbook for Diabetologists, Surgeons, and Diabetes Nurses, as well as the instructions for use of the Accu-Chek DiaPort infusion set and the Accu-Chek DiaPort fixation disc.

**Disposal**  
This set contains pointed or sharp-edged components. Dispose of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured. For more information, contact your local council or authority.

**Reporting of Serious Incidents**  
For a patient/user/third party in the European Union and in countries with identical regulatory regimes, if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

**Summary of Safety and Clinical Performance**  
You find the SSCP after the launch of the European database on medical devices (Eudamed) here: <https://ec.europa.eu/tools/eudamed>

**Distributed in the United Kingdom by:**  
Roche Diabetes Care Limited 3, Roche Avenue, Burgess Hill West Sussex, RH15 9RY, United Kingdom  
Accu-Chek Pump Careline <sup>1)</sup>: UK Freephone number: 0800 731 22 91  
ROI Freephone number: 1 800 68 23 51

**Australia**  
Roche Diabetes Care Australia Pty. Limited  
Insulin Pump Support: 1 800 933 457  
Australia.insulinpumps@roche.com  
www.accu-check.com.au

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Replace the infusion set and the fixation disc according to the intervals specified in the respective instructions for use.

**WARNING**  
**Risk of hyperglycaemia**  
The use of system components beyond the specified replacement intervals can cause occlusions or leaks that can lead to insufficient insulin delivery. Replace system components according to the specified intervals.

**WARNING**  
**Risk of inappropriate therapy**  
The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

**WARNING**  
**Risk of infection**  
A damaged packaging may no longer be sterile or may contain a damaged component. Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

**WARNING**  
**Risk of injury**  
If you modify the system or if you use tools that were not provided with the set, you may injure internal organs and tissues.

**WARNING**  
**Risk of hyperglycaemia and infection**  
The use of sterile components after their Use by date can lead to infections

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