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
ACCU-CHEK®

Made in France







ACCU-CHEK, ACCU-CHEK SPIRIT COMBO and DIAPORT are trademarks of Roche.

© 2020 Roche Diabetes Care

Roche Diabetes Care GmbH
 Sandhofer Strasse 116
 68305 Mannheim, Germany
www.accu-chek.com

00948011054/G Last update: 2020-08 T1002-D

(EN) Accu-Chek DiaPort membrane replacement set

Intended Use

The Accu-Chek DiaPort system is intended for continuous intraperitoneal insulin infusion (CIPI). For this purpose, it shall be used together with an external Accu-Chek Spirit Combo insulin pump. The Accu-Chek DiaPort system delivers insulin from the Accu-Chek insulin pump through the infusion set and the catheter into the peritoneal cavity.

⚠ WARNING
Risk of hyperglycaemia
The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type I or type II diabetes who meet the indication for CIPI therapy. The Accu-Chek DiaPort system is only for use by patients and healthcare professionals specifically trained in the use of the individual components relative to their areas of expertise.

⚠ WARNING
Risk of inappropriate therapy
The Accu-Chek DiaPort system is designed for intended users only. Do not use the system with patients that do not meet the indications, or where one or more contraindications apply.

The Accu-Chek DiaPort membrane replacement set enables the change of the Accu-Chek DiaPort membrane. The membrane must be changed every 6 months to prevent leakages.

⚠ PRECAUTION
Risk of hyperglycaemia
The use of system components beyond the specified replacement intervals can cause occlusions or leaks that can lead to insufficient insulin delivery. Replace system components according to the specified intervals.

The product is intended only for use by healthcare professionals trained in the Accu-Chek DiaPort system, using the tools provided with the set. Membrane replacement must be performed in a clean environment under sterile conditions.

⚠ WARNING
Risk of different types of harm
Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explanation of the system. Do not modify the system and adhere to the instructions for use.

Indications and Contraindications
See the instructions for use of the Accu-Chek DiaPort implantation set.

Component Overview
This is a sterile, single use product that has been sterilised with ethylene oxide.

- A Fork wrench
- B Screw-off screwdriver
- C Screw-in screwdriver with pre-mounted membrane
- D Fixation disc
- E Disc applicator
- F Sticker

⚠ WARNING
Risk of infection
A damaged packaging may no longer be sterile or may contain a damaged component. Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

Risk of suffocation
This product contains small parts. Keep small parts away from people who might swallow them, for example, small children.

Storage
Keep the product dry and away from sunlight. Store the product between 5 and 45 °C.

⚠ WARNING
Risk of hyperglycaemia and infection
The use of sterile components after their *Use by date* can lead to infections and to abscesses. Do not use components that have passed their *Use by date*.

Before You Get Started

Do not use this product for any medication other than U100 short-acting regular human insulin or sterile physiological saline solution. It may be advisable to take a 30 min warm bath prior to the procedure to dissolve leftovers or crystallisation from insulin. When taking off the insulin pump, do not put it in STOP mode. The insulin dripping from the infusion set keeps the tip clean. Read the Accu-Chek DiaPort Handbook for Diabetologists, Surgeons, and Diabetes Nurses, as well as the instructions for use of the Accu-Chek DiaPort infusion set and the Accu-Chek DiaPort fixation disc.

Replacing the Membrane

1. Use the disc applicator to remove the fixation disc. Disinfect the port and the surrounding area using a non-sticky disinfectant. Cover the surrounding area with a sterile surgical drape.
2. Use the fork wrench to hold the port securely.
3. Insert the screw-off screwdriver at a 90° angle into the membrane. Carefully unscrew the membrane by turning the screw-off screwdriver in a counterclockwise direction. When unscrewing the membrane, make sure the screw-off screwdriver remains at a 90° angle. Use the fork wrench to keep the port in place. Take care not to displace the catheter hub from the indent inside the port body.

⚠ WARNING

Risk of hyperglycaemia
A dislodged catheter hub can lead to insufficient insulin delivery. Make sure that the catheter hub is correctly aligned during and after the change of the membrane. The notch of the catheter hub fits into the port body in a defined direction only.

4. Use **only the screw-in screwdriver** with a new, pre-mounted membrane to screw the membrane into the port body until it is hand-tight. Screw-in the membrane in one continuous motion. A continuous motion avoids the detachment or loosening of the pre-mounted membrane in the screw-in screwdriver.
5. Use the disc applicator to position a new fixation disc on the port. Stick the membrane replacement set sticker on the implant card carried by the patient.

⚠ PRECAUTION

Risk of hyperglycaemia
Handling errors or defective components can lead to insufficient insulin delivery. After replacing a component of the system, check the blood glucose level at least once within 1 to 3 hours.

Disposal
This set contains pointed or sharp-edged components. Disposal of pointed or sharp-edged components according to local regulations, in such a way that no third party can be injured. For more information, contact your local council or authority.

Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime; if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

Summary of Safety and Clinical Performance
You find the SSP after the launch of the European database on medical devices (Eudamed) here: <https://ec.europa.eu/tools/eudamed>

Distributed in the United Kingdom by:
Roche Diabetes Care Limited
Charles Avenue, Burgess Hill
West Sussex, RH15 9RY, United Kingdom
Accu-Chek Pump Careline 1:
UK Freephone number: 0800 731 22 91
ROI Freephone number: 1 800 88 23 51

¹⁾ calls may be recorded for training purposes
Some mobile operators may charge for calls to these numbers.
burgesshill.insulinpumps@roche.com
www.accu-chek.co.uk
www.accu-chek.ie

Australia
Roche Diabetes Care Australia Pty. Limited
Insulin Pump Support: 1800 633 457
australia.insulinpumps@roche.com
www.accu-chek.com.au

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Roche Diabetes Care France SAS