# Handbook for Port Users

The Accu-Chek DiaPort system

EN DE FR IT NL NO

<b>EN</b> Handbook for Port Users .....	2-39
<b>DE</b> Handbuch für Portträger.....	40-77
<b>FR</b> Guide pour les utilisateurs du port d'accès péritonéal.....	78-115
<b>IT</b> Manuale per i pazienti con sistema di accesso per l'infusione continua di insulina intraperitoneale .....	116-153
<b>NL</b> Handboek voor infusiepoortgebruikers .....	154-191
<b>NO</b> Håndbok for portbrukere .....	192-229

**Distributed in the United Kingdom by:**

Roche Diabetes Care Limited  
Charles Avenue, Burgess Hill  
West Sussex, RH15 9RY, **United Kingdom**  
Accu-Chek Pump Careline <sup>1)</sup>:  
UK Freephone number: 0800 731 22 91  
ROI Freephone number: 1 800 88 23 51

<sup>1)</sup> calls may be recorded for training purposes  
Some mobile operators may charge for calls to these numbers.

burgesshill.insulinpumps@roche.com  
www.accu-chek.co.uk  
www.accu-chek.ie

**Australia**

Roche Diabetes Care Australia Pty. Limited  
Insulin Pump Support: 1800 633 457  
australia.insulinpumps@roche.com  
www.accu-chek.com.au

# Product Overview

Components of the Accu-Chek DiaPort system:

<p>1</p> 	<p>Handling aid</p>
<p>2</p> 	<p>Infusion set with ball cannula</p>
<p>3</p> 	<p>Disc applicator</p>
<p>4</p> 	<p>Fixation disc</p>

5

Membrane


Port body

6

Catheter

7

Implant card

### Healthcare Professional

Name and address:

---

Phone:

---

Email:

---

### Surgeon

Name and hospital address:

---

Phone:

---

Email:

---

# About this Handbook

This handbook provides you with information and guidance about the Accu-Chek DiaPort system and continuous intraperitoneal insulin infusion (CIPII). The Accu-Chek DiaPort is referred to as “port” when talking about the port body that contains the membrane and as “system” when talking about the port and its accessories, such as, the catheter or the fixation disc.

Together with your healthcare professional, it helps you to:

- Decide to switch to the Accu-Chek DiaPort system or remain with your current therapy.
- Understand what is involved in the implantation of the system, postoperative care and the long term care for the system.
- Understand how managing your diabetes will change.
- Get an idea of daily life with the system after the implantation.

The Accu-Chek DiaPort system is only for use by healthcare professionals and patients specifically trained in the use of the individual components relative to their areas of expertise. The system may only be implanted by a trained surgeon working with a certified Accu-Chek DiaPort center of excellence.

Prior to using this product or any of its components, consult the instructions for use that come with the product. For explanations of the medical terms used in this handbook, see the Glossary section.


# Table of Contents

Product Overview .....	3
About this Handbook.....	5
Introduction .....	7
1 Deciding on the Accu-Chek DiaPort System.....	8
1.1 Indications.....	8
1.2 Intended Users .....	9
1.3 Contraindications .....	9
1.4 General Warnings.....	10
1.5 Advantages and Disadvantages .....	11
1.6 Switching from Subcutaneous Pump Therapy.....	12
2 System Components.....	13
2.1 Port Body.....	13
2.2 Membrane.....	14
2.3 Catheter.....	15
2.4 Fixation Disc .....	16
2.5 Infusion Set.....	17
3 Implantation.....	18
3.1 Surgical Procedure .....	18
3.2 Implantation Site.....	20
3.3 What Will the Implanted Port Look Like?.....	21
3.4 Postoperative Care.....	22
4 Therapy.....	24
4.1 Using an Accu-Chek Spirit Combo Insulin Pump with the System.....	24
4.2 How Does the Fixation Disc Work? .....	26
4.3 How Does the Infusion Set Work? .....	27
5 Daily Life .....	28
5.1 Professional and Leisure Activities.....	28
5.2 Medical Issues.....	33
5.3 Hypoglycaemia and Hyperglycaemia.....	34
5.4 Daily Care.....	35
Glossary .....	37
Symbol Descriptions .....	39

# Introduction

The Accu-Chek DiaPort system enables continuous intraperitoneal insulin infusion: That is, infusion of insulin into the peritoneal cavity with a compatible Accu-Chek insulin pump and a specially designed ball cannula infusion set. The port has a titanium body with a catheter that is placed in the abdomen. The top of the port rises above the surface of the skin by about 5 mm (0.19 inches), while a flower-shaped plate is surgically placed under the skin in a subcutaneous pocket, providing stability for the system.

The system has been designed for people with diabetes who cannot fully benefit from continuous subcutaneous insulin infusion (CSII). When considering continuous intraperitoneal insulin infusion (CIPII) therapy, discuss all aspects with your healthcare professional.

*The Accu-Chek DiaPort system with an Accu-Chek Spirit Combo insulin pump.*

# 1 Deciding on the Accu-Chek DiaPort System

Together with your healthcare professional and an expert from a center of excellence, you are considering the system as a solution for your specific situation. Discuss with him or her the indications and contraindications, as listed below.

You may already be familiar with subcutaneous insulin pump therapy (often referred to as CSII therapy). When you change to the port therapy, the insulin will no longer be infused under your skin but into your abdomen (into the so called peritoneal cavity). The new therapy is called continuous intraperitoneal insulin infusion (CIPII).

## 1.1 Indications

The Accu-Chek DiaPort system may be a suitable solution if one or more of the following indications exist:<sup>1</sup>

- People with diabetes mellitus requiring insulin in whom subcutaneous insulin therapy fails as shown by not achieving glycaemic targets, or only achieving glycaemic targets at the expense of increased rates of severe hypoglycaemia in spite of the use of advanced technology, for instance in subcutaneous insulin resistance.
- People with diabetes mellitus requiring insulin in whom subcutaneous insulin is not feasible because of skin disorders, for instance in the presence of lipohypertrophy or -atrophy.

<sup>1</sup> Schade D.S., Duckworth W.C., *In Search of the Subcutaneous-Insulin-Resistance Syndrome*. N Engl J Med, 1986 Jul 17;315(3):147-53.  
Soudan B., et al., *Extreme Subcutaneous Insulin Resistance: A Misunderstood Syndrome*. Diabetes Metab, 2003 Nov;29(5):539-46.  
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.  
Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.  
Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diaport™-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.  
Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipohypertrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.  
Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.  
Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821.

## 1.2 Intended Users

Intended users of the Accu-Chek DiaPort system are insulin-dependent people with type I or type II diabetes who meet the indication for CIPII therapy.

The user or his/her permanent caregiver must be fully aware of the details of the Accu-Chek insulin pump therapy and he or she must be well-trained in the use of the insulin pump therapy as well as the backup systems.

The user's healthcare professional decides whether a user, either alone or with mandatory and permanent support of a caregiver is able to use an insulin pump for intraperitoneal infusion. Special attention has to be paid when using intraperitoneal insulin pump therapy in children or neonates. There is no experience with intraperitoneal insulin pump therapy in neonates.

### **WARNING**

#### **Risk of inappropriate therapy**

The Accu-Chek DiaPort system is designed for intended users only.

Do not use the system with patients that do not meet the indications, or where one or more contraindications apply.

## 1.3 Contraindications

The Accu-Chek DiaPort system may not be suitable if one or more of the following contraindications apply:

- Patients not considered suitable for using an external insulin pump, or patients with contraindications for CSII:
  - Poor compliance with current treatment
  - Recurrent failure to attend clinics
  - Inability to follow a treatment plan
  - Unwillingness to calculate meal dosages or to perform the necessary number of blood glucose tests per day for this type of therapy
  - Evidence of psychiatric conditions
- Allergy to one or more of the port materials (titanium, polyethylene, etc.) including its infusion set
- Disorders of the gastrointestinal tract interfering severely with the intraperitoneal route of administration (for example, colonic disease)
- Continuous ambulatory peritoneal dialysis (CAPD)
- Concomitant drugs administered by intraperitoneal infusion (for example, chemotherapy)
- Severely impaired immune response
- High levels of anti-insulin antibodies

## 1.4 General Warnings

Do not use this product for any medication other than U100 short-acting regular human insulin or sterile physiological saline solution.

Replace the infusion set and the fixation disc according to the intervals specified in the respective instructions for use.

### **PRECAUTION**

#### **Risk of hyperglycaemia**

If you do not know your current blood glucose level, you may not be able to react in time to some situations, for example, if an occlusion occurs.

Check your blood glucose level at least 4 times a day and whenever you are in doubt about your current blood glucose level.

#### **Risk of hyperglycaemia**

Handling errors or defective components can lead to insufficient insulin delivery.

After replacing a component of the system, check your blood glucose level at least once within 1 to 3 hours.

#### **Risk of hyperglycaemia**

An undetected occlusion can lead to insufficient insulin delivery.

If you experience an unexpected rise in your blood glucose level, check the insulin pump and infusion set for occlusions. If you are not sure whether the infusion set is working properly, replace the infusion set. If you cannot resolve the occlusion, contact your healthcare professional **immediately**.

#### **Risk of hyperglycaemia**

The use of system components beyond the specified replacement intervals can cause occlusions or leaks that can lead to insufficient insulin delivery.

Replace system components according to the specified intervals.

### **WARNING**

#### **Risk of interruption of therapy and scarring**

If the port body sinks into the subcutaneous tissue, explantation may be required.

If you notice that the port body starts sinking, contact your healthcare professional immediately.

#### **Risk of infection**

A damaged packaging may no longer be sterile or may contain a damaged component.

Before using a new component, make sure its packaging is not damaged. If the packaging is damaged, do not use its contents.

#### **Risk of hyperglycaemia and infection**

The use of sterile components after their *Use by* date can lead to infections and to abscesses.

Do not use components that have passed their *Use by* date.

#### **Risk of suffocation**

This product contains small parts.

Keep small parts away from people who might swallow them, for example, small children.

**Risk of different types of harm**

Modification of system components or disregard of the instructions for use may prevent the system from working as intended. This can cause different types of harm, for example, infection, hyperglycaemia, damage to the skin, or scarring due to explantation of the system. Do not modify the system and adhere to the instructions for use.

## 1.5 Advantages and Disadvantages

### Advantages

Therapy using the Accu-Chek DiaPort system can have the following advantages:

- In contrast to subcutaneous injections or subcutaneous pump therapy, the insulin flows through the catheter into the peritoneal cavity and reaches the liver more quickly. The insulin starts taking effect within a few minutes.
- Using this method may improve your diabetes control, especially reducing frequent hypoglycaemia, as CIPII more closely mimics the normal physiological secretion of insulin. CIPII can lead to an improved blood glucose profile.
- Compared to subcutaneous insulin infusion, HbA<sub>1c</sub> values can be improved.<sup>2</sup>
- The frequency of unexpected and rapid drops in the blood glucose level can be reduced.
- You no longer need to change infusion sites every 2 to 3 days, as the insulin is infused directly through the port catheter.

### Disadvantages

Therapy using the Accu-Chek DiaPort system can have the following disadvantages:


- There is a risk of infection of the tissue around the top of the port body, especially just after the implantation. However, taking the recommended precautions after implantation helps to reduce the risk of infections.
- There can be adhesions to the abdominal wall, which may result in pain, or the tip of the catheter may become encapsulated. In this case, it is sometimes necessary to replace the catheter.
- The catheter can become blocked, which is called occlusion.
- In rare cases, the port must be explanted, for example, in case of serious or recurrent non-treatable infections. However, once the infection is cleared, a new port can be implanted at a new site.

<sup>2</sup> Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipodystrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.

Liebl A., et al., *Long-term clinical evaluation of the new Accu-Chek® diapor, a port system for continuous intraperitoneal insulin infusion: 24-month results*. Diabetes, 2014; 63: A241.

**The incidence of severe hypoglycaemia in participants with CIPII was less than half of the incidence with CSII.<sup>3</sup>**

## 1.6 Switching from Subcutaneous Pump Therapy

The switch from subcutaneous insulin pump therapy (CSII) to intraperitoneal insulin therapy (CIPII) can take place immediately after the implantation of the port. However, your healthcare professional decides when your insulin pump can be connected. When the first insulin dose is delivered through CIPII, a healthcare professional must be present.

For certain indications (like subcutaneous insulin resistance), the insulin requirement may drop significantly when switching to port therapy. Your insulin dose is adjusted by your healthcare professional.

With intraperitoneal insulin delivery, there is less delay between insulin delivery and its absorption. The reduced delay affects the calculation and infusion of boluses. If you use a bolus calculator, the parameters “offset time” and “acting time” might need to be adjusted. Discuss these factors with your healthcare professional before starting CIPII therapy.

### **WARNING**

#### **Risk of hyperglycaemia or hypoglycaemia**

When starting CIPII therapy, your insulin needs are likely to change.

Speak to your healthcare professional about the adjustments to your therapy parameters.

<sup>3</sup> Liebl A., et al., *A Reduction in Severe Hypoglycaemia in Type 1 Diabetes in a Randomized Crossover Study of Continuous Intraperitoneal Compared With Subcutaneous Insulin Infusion*. *Diabetes Obes Metab*, 2009 Nov;11(11):1001-8.

## 2 System Components

### 2.1 Port Body

Membrane (1), port body (2), and catheter (3)

The port enables intraperitoneal insulin infusion: That is, infusion of insulin into the peritoneal cavity with an insulin pump and an infusion set. The port has a metal body and a catheter that is placed in the abdomen. The opening of the port rises above the surface of the skin by about 5 mm (0.19 inches), while a flower-shaped plate with a diameter of 26 mm (1.02 inches) is placed under the skin, providing stability for the system.

The system has been re-designed and re-developed with major improvements:

- To make sure that the skin grows correctly around the body of the port, a special felt material is wrapped around the base of the port, acting as an infection barrier.
- The catheter is made of a new material that is more flexible and has a softer and wider tip, to allow free movement in the abdominal cavity. This helps to reduce the risk of adhesions or overgrowing/encapsulation of tissue at the catheter tip.
- The risk of occlusion is expected to be significantly reduced due to an enlarged inner diameter of the catheter and the trumpet-shaped soft tip.
- All the materials used in the port are biocompatible and hypo-allergenic, to reduce the risk of rejection to the implanted components by your body.

#### **PRECAUTION**

##### **Risk of adverse reaction**

In rare situations, the implanted system may cause adverse reactions.

If you experience an adverse reaction, contact your healthcare professional immediately.

#### **WARNING**

##### **Risk of scarring**

In rare situations, the port body may be rejected by the human body, which requires explantation. Scars remain after explantation of the port system.

Materials of the port body:

Component	Materials	Volume	Ratio
Base	Titane 6 AL 4V ELI + Anodisation	630.43 mm <sup>3</sup> (0.038 in <sup>3</sup> )	75.20 %
Glue	Silicone adhesive MED-1511	Traces	Traces
Primer	MED6-161	Traces	Traces
Polyester felt	Polyester felt 6077 (Dacron)	207.94 mm <sup>3</sup> (0.013 in <sup>3</sup> )	24.80 %

Expected lifetime: 5 years

It is recommended not to explant the port after 5 years if the port is still functional.

## 2.2 Membrane

The membrane seals the access to the peritoneal cavity. When the ball cannula is removed, the port is automatically sealed tight by the silicone membrane. No water or solution can go through the port.

The membrane must be changed every 6 months to prevent leakage. This is a virtually pain-free procedure, but must be performed at a center of excellence by a trained healthcare professional or surgeon under sterile conditions.

Materials of the membrane:

Component	Materials	Volume	Ratio
Lid	Titane 6 AL 4V ELI + Anodisation	119.89 mm <sup>3</sup> (0.007 in <sup>3</sup> )	62.67 %
O-ring	SILOPREN LSR 4040 (40 shore A)	9.19 mm <sup>3</sup> (0.00056 in <sup>3</sup> )	4.80 %
Valve	SILOPREN LSR 4060 (60 shore A)	62.02 mm <sup>3</sup> (0.004 in <sup>3</sup> )	32.50 %

Replacement interval: At least every 6 months

## 2.3 Catheter

The catheter enables the infusion of insulin from the external insulin pump into the peritoneal cavity. The catheter is implanted into the abdomen and not visible from the outside of the body. The catheter is flexible, and has a trumpet-shaped soft tip so that it can move around in the abdominal cavity. This free movement helps to reduce the risk of adhesions or overgrowing of tissue on the catheter tip.

3 different lengths are available, 9 cm, 15 cm and 25 cm (3.5 inches, 5.9 inches and 9.8 inches). Your surgeon decides which length is suitable for you.

The catheter must be replaced if it is blocked or overgrown. The replacement must be performed at a center of excellence by a trained surgeon under sterile conditions.

Materials of the catheter:

Component	Materials	Volume	Ratio
Catheter	Tecothane TT2074A-B20	238.05 mm <sup>3</sup> (0.0145 in <sup>3</sup> )	28.30 %
Catheter	Tecothane TT1074A	521.28 mm <sup>3</sup> (0.0318 in <sup>3</sup> )	61.90 %
Catheter	PTFE F polyflon	30.28 mm <sup>3</sup> (0.0184 in <sup>3</sup> )	3.60 %
Hub	Tecothane TT1065D	52.71 mm <sup>3</sup> (0.0032 in <sup>3</sup> )	6.30 %

Expected lifetime: 5 years

It is recommended not to explant the catheter after 5 years if the catheter is still functional.

## 2.4 Fixation Disc

*Fixation disc (1) and disc applicator (2)*

The fixation disc provides stability and flexibility. This is important because the port is not sutured into the subcutaneous tissue. Instead, the port takes hold through progressive in-growth/granulation of the subcutaneous tissue.

*The implanted system with the fixation disc*

Even when the healing process is complete, the fixation disc continues to stabilise the port, absorb impacts, and is also helpful in holding the port when connecting and disconnecting the infusion set. This helps to reduce stress on the tissue surrounding the port site by limiting movement within the implantation site/pocket. Additionally, it helps to prevent the port from sinking or being overgrown by skin. The disc applicator is used for fitting and removing the fixation disc.

You must wear the fixation disc at all times. Change the fixation disc when changing the cartridge of the insulin pump or at least every 6 days. Clean the port and the surrounding tissue with sterile sodium chloride solution or with clean water on a daily basis or at least when changing the fixation disc. Do not use solutions containing alcohol.

### **WARNING**

#### **Risk of interruption of therapy and scarring**

The fixation disc prevents the port body from sinking into the subcutaneous tissue, which may require explantation.

Wear the fixation disc at all times.

Materials of the fixation disc and disc applicator:

Component	Materials	Volume	Ratio
Fixation disc	PP BORMED RF825 MO	n/a	n/a
Disc applicator	Compound based on PP 1013 H1	n/a	n/a
Disc applicator	Blue colorant Lifocolor PP (0.5 %)	n/a	n/a

Replacement interval: At least every 6 days

## 2.5 Infusion Set

### *Infusion set*

The infusion set connects the system with an external Accu-Chek Spirit Combo insulin pump. When using the infusion set for the first time, a healthcare professional must be present.

The infusion set has been specially designed for use with the port, and has a ball cannula that is inserted into the port body membrane, which enables the connection with the catheter. This solution avoids needle-stick injuries, and is virtually pain free. After disconnection of the infusion set, the membrane automatically seals the access to the peritoneal cavity.

*The infusion set with the handling aid attached to the insulin pump*

Replacement interval: At least every 6 days

# 3 Implantation

## 3.1 Surgical Procedure

The port is implanted under general anesthesia during a minor operation that takes approximately 20 to 30 minutes. Discuss the placement of the port, the implantation, and the anesthesia in advance with your healthcare professional, surgeon, and anesthetist.

The operation takes place in a hospital that is associated with the center of excellence. It is recommended that you remain in hospital for a couple of days.

The implantation starts with a small incision approximately 4 cm (1.57 inches) wide and by cutting a small horizontal pocket under the skin. Then, a circular hole for the top outer part of the port (including the membrane) is cut through the skin. The port is placed under the skin with the top part protruding through the prepared hole. The catheter is buried in the peritoneal cavity through a small cut in the abdominal wall inside the pocket.


There may be slight to moderate pain following the operation, but this can be managed with analgesics. If the level of pain seems excessive, or is causing you discomfort, consult your healthcare professional immediately.

### **WARNING**

#### **Risk of abdominal pain**

The catheter of the Accu-Chek DiaPort system is placed in the abdominal cavity. Although the catheter is made of a flexible material, it can cause abdominal pain.

*The implanted Accu-Chek DiaPort system*

## 3.2 Implantation Site

The port can be implanted wherever you prefer in the lower abdominal region as shown by the green shaded area in the pictures, with a few exceptions. It is very important that you take time to find the best place for implanting the port. Your daily activities must not be affected by the positioning of the port. So think carefully about the different activities you undertake in your studies, work or leisure time, especially sporting activities.

Do not position the port on or near pressure points.



For example, the port must **not** be positioned in:

- the region close to the navel.
- the region crossed by a car safety belt.
- the region under or near the trouser belt.
- the region around the rib-arch.
- a place where it can bump against the edge of a table (when seated).

Right-handed patients may prefer to have the port on the left side, and vice versa.

Make sure you discuss a 2<sup>nd</sup> surgical placement option in case there is an issue during surgery, preventing the port from being placed at the original preferred site.

It is recommendable to choose an implantation site in such a way that you can always see the port during handling.

### 3.3 What Will the Implanted Port Look Like?

The operation leaves a scar which is about 4 cm (1.57 inches) long. The exact size depends on the surgical requirements during the operation (for example, unexpected complications).

The top of the port (diameter of 9 mm (0.35 inches)) protrudes above the surface of the skin by up to 5 mm (0.19 inches).

The picture shows the port with the fixation disc and the infusion set several weeks after surgery.

*The implanted Accu-Chek DiaPort system*

## 3.4 Postoperative Care

Normally, the port can be used immediately following implantation, but this decision is under the responsibility of your healthcare professional.

The scar left by the operation, needs about 10 days to heal. Healing of the tissue surrounding port implantation site normally takes between 4 to 6 weeks, if no infection or no other problems occur.

During the first few days after implantation, it is recommended to rest as much as possible. This supports the healing of the wound.

It is very important to keep the port area clean and not to undertake activities that affect the healing of the wound. It is strongly recommended not to lift or carry heavy objects during the healing time. Before you return to your normal daily activities, consult your healthcare professional to ensure it is alright to do so.

### **Dressing the port site and the surgical wound**

The success of the port therapy depends on how well the tissue grows around the port, and how clean the port is kept.

If possible, do not remove the fixation disc for the first couple of days. Wash and dry the fixation disc daily. Replace the fixation disc with a new one every 6 days. It is recommended not to take a shower at least during the first 3 days after implantation. During the healing period, clean the port and the surgical wound with sterile water or sterile physiological saline solution. Then, dress the port with a dry, sterile dressing (padded if necessary) once a day until the wound has healed completely.

After the healing period, it is recommended to continue to take care of the port on a daily basis. If possible, cover the port with a sterile adhesive dressing (not water tight) as this helps to keep the port as clean as possible and reduces the risk of infections.

Take care when washing the area of the wound and the port implantation site, do not rub or apply excessive force. Make sure to clean and dry the fixation disc and the skin. Do not use solutions containing alcohol, as these are too aggressive for the skin, and may slow down the healing process.

Do not graze or pick the scab growing on the wound around the port. This is part of the usual healing process.

## Precautions

- Wash and dry your hands thoroughly before and after touching any components of the port.
- Always observe strict rules of hygiene and cleanliness around the port. Always wipe from clean to dirty areas to prevent the risk of infections.
- Avoid exerting pressure on the port: For example, from belts, tight clothing, or car safety belts.
- After showering or bathing, always pat the skin dry around the port thoroughly.
- Do not place cream around the port area.
- Do not expose the port to direct sunlight or high temperatures for any length of time, as there is a risk of skin burns from the metal components.

### **WARNING**

#### **Risk of skin burns**

If the port body is exposed to high temperatures, for example, in a sauna, skin burns can occur. Do not expose the port body to high temperatures.

# 4 Therapy

## 4.1 Using an Accu-Chek Spirit Combo Insulin Pump with the System

The switch from CSII to CIPII normally takes place immediately after the implantation of the port. However, your healthcare professional will advise you when he or she plans to switch. Depending on your indication for CIPII, your daily insulin need may drop significantly when switching to port therapy. This point also needs to be discussed with your healthcare professional.

Using an Accu-Chek Spirit Combo insulin pump with the port is very similar to using an insulin pump for continuous subcutaneous insulin infusion. Not all types of insulin are suitable for CIPII therapy. **Use only U100 short-acting regular human insulin.** Discuss this with your center of excellence team. We recommend using **Insuman® Infusat** because it showed a low tendency towards occlusions after several weeks.<sup>4</sup>

### **WARNING**

#### **Risk of hyperglycaemia**

The use of a non-compatible insulin pump can lead to insufficient insulin delivery. Use the Accu-Chek DiaPort system only with an Accu-Chek Spirit Combo insulin pump.

The main difference with port therapy is the insulin action: The insulin acts more rapidly and effectively than with subcutaneous therapy.<sup>5</sup> The faster action of insulin enabled by the port means that, for example, you can administer pre-meal boluses just before you sit down to eat. When the meal has a high fat content, it is a good option to use an extended bolus or a multiwave bolus. This means that together with your healthcare professional you will have to reprogram the parameters of the bolus calculator. The settings must be adapted accordingly. Discuss the use of these settings with your healthcare professional.

Stopping the pump immediately affects the availability of insulin. This gives you a high level of control of your blood glucose levels and can help to reduce the risk of hypoglycaemic events.

<sup>4</sup> Liebl A., et al., *Evaluation of the New ACCU-CHEK® DIAPORT, a Port System for Continuous Intraperitoneal Insulin Infusion, in Patients With Type 1 Diabetes: First 6-Month Results.* Diabetes, 2013; 62: A247-A248.

<sup>5</sup> Schade D.S., et al., *Normalization of plasma insulin profiles with intraperitoneal insulin infusion in diabetic man.* Diabetologia 1980; 19(1):35-39. Radziuk J., et al., *Splanchnic and Systemic Absorption of Intraperitoneal Insulin Using a New Double-Tracer Method.* Am J Physiol, 1994 May;266 (5 Pt 1):E750-9.

Selam J.L., *Implantable Insulin Pumps.* Lancet, 1999 Jul 17;354(9174):178-9.

Catargi B., et al., *Comparison of Blood Glucose Stability and HbA1C Between Implantable Insulin Pumps Using U400 HOE 21PH Insulin and External Pumps Using Lispro in Type 1 Diabetic Patients: A Pilot Study.* Diabetes Metab, 2002 Apr;28(2):133-7.

## How long can the insulin pump remain disconnected?

Compared to CSII therapy, no subcutaneous insulin depot is present. This means that a few minutes after the insulin pump has been disconnected, very little insulin is still available. Discuss this point with your healthcare professional.

### **WARNING**

#### **Risk of hyperglycaemia**

Compared to CSII therapy, CIPII therapy does not have a subcutaneous insulin depot. If you remove the insulin pump, very little insulin remains in your body, which may increase the risk of hyperglycaemia.

Ask your healthcare professional for how long you can remove the insulin pump safely.

When you reconnect the insulin pump to the port, check your blood glucose level and if necessary, administer a bolus.

If you disconnect your pump for a short time, do **not** put it in STOP mode, for example, when taking a shower. The insulin dripping from the infusion set keeps the tip clean. Use the handling aid to attach the infusion set on the luer-lock connector of the pump. Do not touch the ball cannula directly after removal of the protective cap.

## What to do if your blood glucose level changes unexpectedly?

First, eliminate all reasons that may be related to the insulin pump or the infusion set.

If an occlusion alarm occurs after changing the infusion set and it is clear that the problem is not related to the insulin pump or to the infusion set, then the catheter is probably blocked.

If you suspect that the catheter is blocked, **immediately** contact your healthcare professional and switch to another insulin delivery method. If the catheter is blocked, a catheter replacement may be required. Replacement of the catheter is performed at the center of excellence.

## What to do if the port or insulin pump are damaged or are not working properly?

Immediately contact your healthcare professional and switch to another insulin delivery method.

## 4.2 How Does the Fixation Disc Work?

The fixation disc is fitted after the implantation of the port and must be worn at all times. The fixation disc holds the port body in place and helps to ensure the initial healing of the wound by providing stability. This is important because the port is not sutured into the subcutaneous tissue, but instead takes hold through progressive in-growth of the surrounding tissue. Even when the healing process is complete, the fixation disc continues to stabilise the port, absorbs impacts, and is also helpful in holding the port steady when connecting and disconnecting the infusion set.

### **WARNING**

#### **Risk of interruption of therapy and scarring**

The fixation disc prevents the port body from sinking into the subcutaneous tissue, which may require explantation.

Wear the fixation disc at all times.

For detailed handling instructions, refer to the instructions for use provided with the product.

## 4.3 How Does the Infusion Set Work?

The infusion set has a ball cannula, which is specially designed for the port. The ball cannula is connected directly to the catheter within the port body by being inserted through the membrane. The ball cannula also avoids the risk of needle-stick injury.

Do not touch the ball cannula after you removed the protective cap. Always use the handling aid to operate the infusion set.

### **PRECAUTION**

#### **Risk of hyperglycaemia**

The infusion set can come loose if pulling forces act on the tubing.

Create a safety loop with the tubing of the infusion set to prevent loosening the infusion set by accidentally pulling on the tubing.

#### **Risk of strangulation**

If the tubing of the infusion set is wrapped around the neck, there is a risk of strangulation.

Pay attention to the tubing of the infusion set.

The fill amount of the infusion set depends on its length:

Infusion set length	70 cm (27.56 inches)	100 cm (39.37 inches)
Fill amount	13 U (U100 insulin)	17 U (U100 insulin)

### **WARNING**

#### **Risk of hypoglycaemia**

When filling the infusion set with insulin, the insulin pump delivers a considerable insulin amount at once.

Do not fill the infusion set while it is connected to the body. You may deliver an uncontrolled insulin amount.

#### **Risk of hyperglycaemia**

If you tighten the luer-lock connector too firmly, it may become damaged, which can lead to insufficient insulin delivery.

Do not overtighten the luer-lock connector.

For detailed handling instructions, refer to the instructions for use provided with the product.

# 5 Daily Life

## 5.1 Professional and Leisure Activities

Your day-to-day activities and therapy with an Accu-Chek Spirit Combo insulin pump should not be affected by the port system. However, if in doubt, consult your healthcare professional.

### **Continuing work**

In the first weeks after implantation, you may need to reduce your normal daily activities until the healing process is completed. The extent to which you may need to reduce your activities is determined by the type of job you do. Except for manual labor or work which may involve risks associated with contact or pressure on the port system, the port should have little impact on your professional activities.

It is particularly important, in the weeks following the operation, that you avoid contact or pressure on the port, so that the wound can heal properly.

### **Driving a vehicle**

When choosing the site for the implantation together with your healthcare professional and the surgeon, take into account that you may have to drive a vehicle. Therefore, it is crucial to position the port so that the safety belt does not interfere with the port or the wearing of the insulin pump. However, avoid exerting pressure on the port from tight clothing or from clothing belts. Before driving, follow the same rules as with CSII therapy. There are no differences in terms of diabetes therapy requirements.

**Physical work**

Once the healing process of the subcutaneous pocket and tissue surrounding the port body is complete, generally after 4 to 6 weeks, there should be no problem with lifting or carrying weights.

**Sport**

Participating in sporting activities requires some care and common sense. The first rule is that the healing process must be complete, with no significant inflammation. In addition, the port must not be subjected to strong pressure or strain.

Activities such as hiking or cycling generally do not present difficulties. However, it is not recommended to participate in contact sports with risk of injury, such as football or boxing.

## Swimming

Swimming is permitted, provided that you only swim in pools with chlorinated water, clean lake water, or in the sea. As soon as the infusion set is removed, the membrane of the port is watertight.

Before going for a swim:

1. Disconnect the infusion set from the port.  
Note that you should discuss with your healthcare professional for how long you can disconnect the infusion set safely. For more information, see chapter 4.
2. It is recommended to cover the port with a waterproof adhesive dressing.  
Make sure the adhesive dressing is large enough so that it does not stick to the membrane or the port body.

After your swim:

1. Remove the waterproof adhesive dressing.
2. Clean the skin surrounding the port with sterile or clean water, and allow it to dry completely.
3. Connect the infusion set to the port.
4. Check your blood glucose level and if necessary, administer a bolus.

## Sunbathing/Solarium

Do not expose the port to direct sunlight for any length of time, as the port body can become very hot and burn the skin around the port. This may affect the stability of the port and can be a cause of infection.

### **WARNING**

#### **Risk of skin burns**

If the port body is exposed to high temperatures, for example, in a sauna, skin burns can occur. Do not expose the port body to high temperatures.

If you want to take a sunbath or go to a solarium, attach a protective dressing over the area of the port. The dressing must be placed in such a way that the port does not heat up to a point where it can burn the surrounding skin.

## Showering

During the 3 days after implantation, it is not recommended to take a shower. Afterwards, when showering, it is recommended that you take this opportunity to clean the port and the surrounding skin.

During the healing period: After your shower, the port and the surgical wound must be dressed with a dry, sterile dressing (padded if necessary) until the wound has healed completely.

After the healing period: After your shower, always dry the port very carefully with a clean towel, without pressure or rubbing.

## Sexual intercourse

There are no additional requirements. Take general care.

## Sauna

Sauna is not recommended, as the port body can become very hot and burn the surrounding skin. Sauna also increases the risk of infection.

## **Sleeping on the abdomen**

Sleeping on the abdomen seems not to be an issue as reported by patients.

## **Flying**

Travel by aircraft does not present problems that are specific to the port. The port is detected by security portals because it is made with titanium, which is already used for many implanted devices.

With the help of your implant card, you can explain what the port is. Normal precautions for the insulin pump still apply. For more information, refer to the User's Manual of the insulin pump.

## **Travelling**

If you have been using continuous subcutaneous insulin infusion (CSII) until now, you can manage travelling in a similar way. However, remember that situations may arise in which you may need specialised medical attention: For example, if the catheter is blocked. Therefore, discuss a backup plan with your healthcare professional before you travel.

If you are not sure about hygienic conditions at your destination, cover your port with an adhesive dressing. Discuss this with your healthcare professional.

In any case, take extra supplies with you, including insulin pens or needles to switch from port therapy in case of a malfunction of the port.

## 5.2 Medical Issues

### Visits to your healthcare professional

As long as the therapy with the port is working well, the same rules as for continuous subcutaneous insulin pump therapy are applied. Discuss the frequency of your visits with your healthcare professional. Because the insulin acts more rapidly and effectively after starting CIPII therapy, immediately contact your healthcare professional whenever complications or difficulties arise.

### Membrane replacement

Visit your healthcare professional at the center of excellence at least every 6 months for changing the membrane. The membrane replacement must be performed by a healthcare professional in a clean environment and under sterile conditions.

Prior to membrane/catheter replacement, it may be advisable to take a 30 minutes warm bath to dissolve leftovers or crystallisation from insulin.

### Pregnancy

Pregnancies with the port are known without complications.<sup>6</sup> No removal of the port was necessary.

If you become pregnant or if you are planning to do so, discuss this with your healthcare professional.

<sup>6</sup> Fonseca V.A., et al., *Diabetic Pregnancy Managed With Intraperitoneal Insulin*. Diabet Med, Jan-Feb 1987;4(1):74-6.

Fonseca V.A., et al., *Successful Pregnancy in Diabetic Controlled With Intraperitoneal Insulin*. Diabetes Care, Jul-Aug 1987;10(4):541-2.

Schnell O., et al., *A case of diabetic pregnancy controlled with a percutaneous access device for intraperitoneal insulin infusion*. Diabetes Care. 1994;17(11):1354-1355.

## 5.3 Hypoglycaemia and Hyperglycaemia

### **Mild hypoglycaemia**

If you experience symptoms of mild hypoglycaemia, proceed in the same way as with CSII or as directed by your healthcare team.

### **Severe hypoglycaemia**

If you experience symptoms of severe hypoglycaemia, the quickest and safest action is to interrupt insulin delivery immediately. You do this by either stopping the insulin pump or removing the ball cannula from the port. Proceed in the same way as with CSII or as directed by your healthcare team. Immediately contact your healthcare professional.

### **Hyperglycaemia and/or diabetic ketoacidosis**

In case of hyperglycaemia, follow the same guidelines as for CSII therapy. First, eliminate all reasons that may be related to the insulin pump or the infusion set. Make sure the infusion set is properly connected and that there are no air bubbles in the tubing. Monitor your blood glucose levels more frequently to see if the corrective action has fixed the issue. If you suspect that the catheter is blocked, immediately contact your healthcare professional and switch to another insulin delivery method.

## 5.4 Daily Care

### Prevention of infection and inflammation

It is important to avoid any infection or inflammation around the port. This is one of the most frequent causes for explantation of the port.

To prevent infection and inflammation of the implantation site, take note of the following points:

- Always handle the port with clean dry hands.
- Check the site of the port once a day, to make sure that everything is alright.
- Do not place cream around the port area.
- Clean the port and the surrounding area with sterile or clean water, or with a physiological saline solution. In general it is not necessary to use a skin disinfectant. Do not use solutions containing alcohol, as these can be too aggressive for the skin. If in doubt, discuss the cleaning of the area with your healthcare professional.
- Carefully dry the skin around the port with a clean towel, without pressure or rubbing.
- Always cover the area with an adhesive dressing if there are risks of soiling or infection. Make sure the adhesive dressing is large enough so that it does not stick to the membrane or the port body.

**Any inflammation or infection must be treated immediately** under supervision of your healthcare professional. As soon as any sign of (even mild) inflammation or infection occurs, immediately inform your healthcare professional and apply a mild disinfectant around the port. A culture will need to be done to identify the pathogenic agent to ensure the correct type of antibiotic therapy is commenced.

#### **WARNING**

##### **Risk of scarring**

In some situations, the port body must be explanted, for example, if medical complications occur. Scars remain after explantation of the port system.

### Care of the infusion set

Avoid contamination of the ball cannula of the infusion set. Use the handling aid to attach the infusion set on the luer-lock connector of the pump. If it does become contaminated, replace the infusion set immediately.

Change the infusion set when changing the cartridge of the insulin pump or at least every 6 days. Change the cartridge of the insulin pump in accordance with the insulin manufacturer's recommendations or at least every 6 days.

For detailed handling instructions, refer to the instructions for use provided with the product.

### Care of the fixation disc

Wash and dry the fixation disc daily. Replace the fixation disc with a new sterile one at least every 6 days.

For detailed handling instructions, refer to the instructions for use provided with the product.

### Emergency situations

If you are admitted to hospital in an emergency (not related to the port) or as a result of an accident, it is important that the hospital staff is aware that you are on CIPII therapy, including the intraperitoneal catheter in your abdomen.

After the implantation, you receive an implant card containing information about the batch number and length of catheter you are using. Always keep the implant card with you.

### Hospitalisation

If you have to be hospitalised, even for non-diabetes therapy related reasons, inform the hospital team about your port therapy.

### Medical Scans and X-Rays

The catheter is visible on the X-Ray or Magnetic Resonance Imaging picture. Take off the insulin pump and remove the insulin pump from the room before the X-Ray examination is performed to avoid any electronic error. See the User's Manual of your insulin pump for warnings and cautions regarding magnetic fields.

Nuclear medicine and ultrasound scans should not affect the port, nor be affected by the port.

# Glossary

<b>Allergy</b>	Over-sensitive reaction to certain materials
<b>CIPII</b>	Continuous intraperitoneal insulin infusion
<b>CSII</b>	Continuous subcutaneous insulin infusion
<b>Catheter</b>	A flexible tube made of synthetic material that delivers the insulin into the abdominal (peritoneal) cavity.
<b>Disc applicator</b>	A tool which is supporting the connection and disconnection of the fixation disc to the body of the port.
<b>Fixation disc</b>	A plastic ring that stabilises the port. It should be kept in place at all times.
<b>Handling aid</b>	An aid that helps to connect or to disconnect the infusion set from the port.
<b>Hyperglycaemia</b>	High blood glucose level
<b>Hypoglycaemia</b>	Low blood glucose level
<b>Infusion set with ball cannula</b>	A special connecting circuit between the insulin pump and the port
<b>Insulin resistance</b>	A sharp drop in the therapeutic effect of insulin, or its failure to work due to insensitivity of the body to insulin
<b>Lipodystrophy</b>	A pathological increase or decrease in the fatty tissue
<b>Percutaneous port system</b>	A port system, where some parts are inside and outside of the body.
<b>Peritoneal</b>	Concerning the space in the abdomen lined with the peritoneum
<b>Physiological effect</b>	An effect that corresponds to the normal functioning of an organism.

---

**Safety loop**

The tubing of the infusion set is taped to the skin near the infusion site in the shape of a loop to avoid dislodgement of the cannula due to accidental dragging on the tubing.

---

**Subcutaneous tissue**

Tissue under the surface of the skin

# Symbol Descriptions

	Consult instructions for use or consult electronic instructions for use
	Caution, refer to safety-related notes in the instructions for use accompanying this product.
	Temperature limit
	Keep away from sunlight
	Keep dry
	Use by (YYYY-MM-DD)
	Do not use if package is damaged and consult instructions for use
	Sterilised using ethylene oxide
	Single sterile barrier system with protective packaging outside
	Double sterile barrier system with protective packaging outside
	Use only once
	Date of manufacture
	Medical device
	Manufacturer
	Unique device identifier
	Catalogue number
	Batch code
	Complies with the provisions of the applicable EU Legislation

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Komponenten des Accu-Chek DiaPort Systems:

<p>1</p>  A blue, dome-shaped plastic component with horizontal ridges and a small protrusion at the bottom. The word "ACCUCHEK" is embossed on the side.	<p>Platzierungshilfe</p>
<p>2</p>  A clear plastic infusion set with a curved tube and a small metal ball cannula at the end.	<p>Infusionsset mit Kugelkanüle</p>
<p>3</p>  A blue, flexible plastic applicator with a curved, ribbed design.	<p>Scheibenappikator</p>
<p>4</p>  A clear plastic circular disc with a central hole and a small protrusion.	<p>Fixationsscheibe</p>

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
# Produktoversikt

Accu-Chek DiaPort-systemets komponenter:

<p>1</p>  A blue, dome-shaped plastic component with a small protrusion at the bottom. The word "ACCUCHEK" is embossed on the side.	<p>Plasseringshjelp</p>
<p>2</p>  A clear plastic infusion set with a curved tube and a small metal ball cannula at the end.	<p>Infusjonssett med kule-kanyle</p>
<p>3</p>  A blue, flexible plastic applicator with a curved, ribbed design.	<p>Skiveapplikator</p>
<p>4</p>  A clear plastic circular disc with a central hole and a small protrusion.	<p>Fikseringsskive</p>

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
Riveline J.P., et al., *Subcutaneous Insulin Resistance Successfully Circumvented on Long Term by Peritoneal Insulin Delivery From an Implantable Pump in Four Diabetic Patients*. Diabetes Metab, 2005 Nov;31(5):496-8.

Baillot-Rudoni S., et al., *Implantable Pump Therapy Restores Metabolic Control and Quality of Life in Type 1 Diabetic Patients With Buschke's Nonsystemic Scleroderma*. Diabetes Care, 2006 Jul;29(7):1710.

Hilgard D., Laengler A., *Severe subcutaneous insulin infusion resistance in a 13 year old child - a case-report of successfully performed therapy with a diapor<sup>TM</sup>-system*. Pediatric Diabetes, 2009; 10(Suppl 11):37.

Noud M.N., et al., *Benefits of intra-peritoneal insulin administration in a child with severe insulin-induced lipodystrophy*. Pediatric Diabetes, 2009;10 (Suppl 11):100.

Lee S.I., Narendran P., *Intraperitoneal Insulin Therapy for a Patient With Type 1 Diabetes With Insulin Injection Site Inflammation*. BMJ Case Reports, 2014 Aug 21.

Garcia-Verdugo R., Erbach M., Schnell O., *A New Optimized Percutaneous Access System for CIPII*. J Diabetes Sci Technol. 2017 Jul;11(4):814-821.

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© 2021 Roche Diabetes Care
ACCU-CHEK, ACCU-CHEK SPIRIT COMBO and DIAPORT are trademarks of Roche.
All other product names and trademarks are the property of their respective owners.
Roche Diabetes Care GmbH