**Rx only**

Made in Switzerland

ACCU-CHEK, ACCU-CHEK FLEXLINK, ACCU-CHEK LINKASSIST and AKKY-CHEK are trademarks of Roche.

© 2019 Roche Diabetes Care

**In US, distributed by:**  
Roche Diabetes Care, Inc.  
Indianapolis, IN 46256, **USA**  
Toll Free 1-800-280-7801

 Roche Diabetes Care GmbH  
Sandhofer Strasse 116  
68305 Mannheim, Germany  
[www.accu-chek.com](http://www.accu-chek.com)

**ACCU-CHEK® LinkAssist**

Insertion device






Instructions for Use

Last update: 2019-09

04818300001/K

**ACCU-CHEK®**


<b>(EN)</b> Instructions for Use	4	<b>(SL)</b> Navodila za uporabo	118
<b>(DE)</b> Gebrauchsanweisung	12	<b>(RO)</b> Instrucțiuni de utilizare	126
<b>(FR)</b> Instructions d'utilisation	20	<b>(RU)</b> Инструкции по использованию	132
<b>(IT)</b> Istruzioni per l'uso	28	<b>(UK)</b> Інструкція із застосування	140
<b>(ES)</b> Instrucciones de uso	34	<b>(ET)</b> Kasutusjuhend	148
<b>(PT)</b> Instruções de utilização	42	<b>(LV)</b> Lietošanas pamācība	154
<b>(NL)</b> Gebruiksaanwijzing	50	<b>(LT)</b> Naudojimo taisyklės	160
<b>(DA)</b> Brugsanvisning	56	<b>(AR)</b> إرشادات الاستخدام	173
<b>(NO)</b> Bruksanvisning	62		
<b>(SV)</b> Bruksanvisning	68		
<b>(FI)</b> Käyttöohje	74		
<b>(EL)</b> Οδηγίες χρήσης	80		
<b>(PL)</b> Instrukcja obsługi	88		
<b>(CS)</b> Návod k použití	94		
<b>(SK)</b> Návod na použitie	100		
<b>(HE)</b> הוראות שימוש	111		
<b>(HU)</b> Használati útmutató	112	<b>(US)</b> Warranty	174

<b>Instructions for Use</b>	<b>4</b>	<b>en</b>
<b>Gebrauchsanweisung</b>	<b>12</b>	<b>de</b>
<b>Instructions d'utilisation</b>	<b>20</b>	<b>fr</b>
<b>Istruzioni per l'uso</b>	<b>28</b>	<b>it</b>
<b>Instrucciones de uso</b>	<b>34</b>	<b>es</b>
<b>Instruções de utilização</b>	<b>42</b>	<b>pt</b>
<b>Gebruiksaanwijzing</b>	<b>50</b>	<b>nl</b>
<b>Brugsanvisning</b>	<b>56</b>	<b>da</b>
<b>Bruksanvisning</b>	<b>62</b>	<b>no</b>
<b>Bruksanvisning</b>	<b>68</b>	<b>sv</b>
<b>Käyttöohje</b>	<b>74</b>	<b>fi</b>
<b>Οδηγίες χρήσης</b>	<b>80</b>	<b>el</b>
<b>Instrukcja obsługi</b>	<b>88</b>	<b>pl</b>
<b>Návod k použití</b>	<b>94</b>	<b>cs</b>
<b>Návod na použitie</b>	<b>100</b>	<b>sk</b>
<b>הוראות שימוש</b>	<b>111</b>	<b>he</b>
<b>Használati útmutató</b>	<b>112</b>	<b>hu</b>
<b>Navodila za uporabo</b>	<b>118</b>	<b>sl</b>

## Overview

**A** Release button

**B** Safety catch

**C** Pretensioning element

**D** Base

**E** Contact areas

## Intended Use

The Accu-Chek LinkAssist insertion device is intended for the insertion of the head set of the Accu-Chek FlexLink infusion set (USA: Accu-Chek Ultraflex infusion set) into the skin. The insertion device can be reused multiple times. The insertion device is intended to be used by one person only and must not be used by other people.

## Before You Get Started

Read these instructions for use and the instructions for use of the Accu-Chek FlexLink infusion set (USA: Accu-Chek Ultraflex infusion set) before using the insertion device.

If the insertion device was dropped or exposed to other mechanical stress, check it for damages, e.g. cracks. Do not use the insertion device if it is damaged.

Information describing a situation for which there is a foreseeable serious hazard with the use of the device, is preceded by the following heading:

### **WARNING**

Information regarding any special care you should exercise for the safe and effective use of the device, or to avoid damage to the device that could occur as a result of use, including misuse, is preceded by the following heading:

### **PRECAUTION**

## Inserting the Cannula into the Skin

### PRECAUTION

#### Risk of injury

If a head set is inserted into the insertion device, the head set could be released accidentally. Do not point the inserted head set at your face or other people.

1

- Push the pretensioning element of the insertion device into the base.

2

- Remove the head set of the infusion set from the packaging.
- Hold the head set by the needle cover.
- Insert the head set, with the blue handling aid pointing upwards, into the insertion device, until you feel the head set click into place.

3

### WARNING

#### Risk of infection or injury

If you remove the needle cover, the introducer needle of head set is not protected anymore. Take care and do not touch the introducer needle.

- Remove the needle cover.

**Tip:** The small tab on the protective film of the adhesive pad indicates the side of the head set to which you will later connect the transfer set. See illustration .

4

- Remove both parts of the protective film from the adhesive pad.

5

- Pull the pretensioning element upwards, away from the base, as far as it will go.

### PRECAUTION

#### Risk of injury

If the safety catch is in the  position, the insertion device is no longer locked. The head set could be released accidentally.

Do not slide the safety catch to the  position, until you have the insertion device placed on your skin.

6

- Place the insertion device onto the chosen infusion site, holding it horizontally. It must be positioned in such a way that the whole surface of each contact area on the base of the insertion device is resting on the skin.

**Tip:** You should be able to comfortably reach the side of the head set which had the small tab on the protective film. That is where you will later connect the transfer set.

- Slide the safety catch to the  position. The insertion device is now unlocked.

7

- Press the release button. The cannula is inserted under the skin.
- Remove the insertion device.

If the head set of the infusion set is still in the insertion device, push the pretensioning element back into the base. Continue with step .

If the cannula is not properly positioned under the skin, remove the adhesive pad from the skin

and pull out the cannula. Start at step **1** with a new head set.

---

**8**

- Press the adhesive pad onto the skin.
- 

**9**

- Continue with removing the handling aid of the head set. Follow the instructions for use of the infusion set.

## **Cleaning the Insertion Device**

- Moisten a cloth with water, with a mild hand soap or with a standard disinfectant (for example, 70 % ethanol).
- Wipe the outside of the insertion device. Do not dip the insertion device into the cleaning fluids.
- Allow the insertion device to dry completely.

## **Storing the Insertion Device**

Keep the product dry and away from sunlight.

Store the insertion device with the pretensioning element retracted and without tension. See illustration **B**.

- Pull the pretensioning element upwards, away from the base, as far as it will go.
- Slide the safety catch to the **I** position.
- Press the release button.

## Disposing of the Insertion Device

- Dispose of the insertion device as plastic waste.

For information about correct disposal, contact your local council or authority.

## Reporting of Serious Incidents

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation 2017/745/EU on medical devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

## Guarantee

There is a guarantee of 18 months on the Accu-Chek LinkAssist insertion device for materials and workmanship, from the date of purchase.

US

***For detailed warranty and service information please refer to the Warranty section at the end of this manual.***

## Customer Support

### **In US, distributed by:**

Roche Diabetes Care, Inc.  
Indianapolis, IN 46256, **USA**  
Toll Free: 1-800-280-7801  
Accu-Chek Customer Care Service Center:  
1-800-688-4578  
[www.accu-chek.com](http://www.accu-chek.com)

### **Australia**

Roche Diabetes Care Australia Pty. Limited  
Pump Support: 1800 633 457  
[www.accu-chek.com.au](http://www.accu-chek.com.au)

### **Singapore**

Accu-Chek ExtraCare line: 6272 9200  
[www.accu-chek.com.sg](http://www.accu-chek.com.sg)

Roche Diabetes Care South Africa (Pty)  
Ltd.

Hertford Office Park, 90 Bekker Road  
Vorna Valley, Midrand

### **South Africa**


1686  
Email: [info@accu-chek.co.za](mailto:info@accu-chek.co.za)  
Call Toll Free: 080-34-22-38-37 (SA only);  
+254 20 523 0560 (Kenya only);  
+27 (11) 504 4677 (Other countries)

### **Distributed in the United Kingdom by:**

Roche Diabetes Care Limited  
Charles Avenue, Burgess Hill  
West Sussex, RH15 9RY, **United Kingdom**  
Accu-Chek Pump Careline <sup>1)</sup>:  
UK Freephone number: 0800 731 22 91  
ROI Freephone number: 1 800 88 23 51  
<sup>1)</sup> calls may be recorded for training purposes  
Some mobile operators may charge for calls  
to these numbers.  
[burgesshill.insulinpumps@roche.com](mailto:burgesshill.insulinpumps@roche.com)  
[www.accu-chek.co.uk](http://www.accu-chek.co.uk)  
[www.accu-chek.ie](http://www.accu-chek.ie)

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## Accu-Chek LinkAssist Limited 18-Month Warranty

Roche Diabetes Care, Inc. ("Roche") warrants to the original purchaser of the Accu-Chek LinkAssist insertion device that your Accu-Chek LinkAssist insertion device will be free from defects in materials and workmanship for 18 months from the date of purchase. If, during this 18-month period, the insertion device does not work properly because of a defect in materials or workmanship, Roche will replace it with a new Accu-Chek LinkAssist insertion device or equivalent product free of charge. The warranty on the replacement insertion device will expire on the date of the original warranty expiration or 90 days after the shipment of a new insertion device, whichever period is longer. The purchaser's exclusive remedy with respect to the Accu-Chek LinkAssist insertion device shall be replacement.

**This warranty does not apply to the performance of an insertion device that has been damaged by accident or has been altered, misused, tampered with, or abused in any way.**

**THE ABOVE WARRANTY IS EXCLUSIVE OF ALL OTHER WARRANTIES, AND ROCHE MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WITHOUT LIMITATION, THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL ROCHE BE LIABLE TO THE PURCHASER OR ANY OTHER PERSON FOR ANY INCIDENTAL, CONSEQUENTIAL, INDIRECT, SPECIAL, OR PUNITIVE DAMAGES ARISING FROM OR IN ANY WAY CONNECTED WITH THE PURCHASE OR OPERATION OF THE ACCU-CHEK LINKASSIST INSERTION DEVICE OR ITS PARTS. NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, IF ANY IS IMPLIED FROM THE SALE OF THE INSERTION DEVICE, SHALL EXTEND FOR A LONGER DURATION THAN EIGHTEEN MONTHS FROM THE DATE OF PURCHASE.**

Some states do not allow limitations on how long an implied warranty will last or the exclusion of incidental or consequential damages, so the above limitation and exclusion may not apply to you. This warranty gives you specific legal rights, which vary from state to state.

## **Warranty and Service Instructions**

All requests for return of the Accu-Chek LinkAssist insertion device under the above warranty must be made to Roche Diabetes Care, Inc. Indianapolis, IN 46256, USA. You will be mailed a return authorization label, which must be affixed to your carton for shipping the system to Roche. Cartons received without this label will be returned to you at your expense.

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